

## 16S Metabarcoding

### Materials:

- 12596-4 Qiagen UltraClean 96 PCR Cleanup Kit (384 samples)
- Q32854 (Thermo Fisher) Invitrogen Qubit dsDNA HS Assay Kit (500 assays)
- MS-102-2002 (Illumina) MiSeq Reagent Kit v2 (300-cycle)
- AM9932 (Thermo Fisher) Ambion Nuclease-free water (not DEPC-treated)
  - pH 8.0, decanted in falcon tube under fume hood
- 13000013 (Thermo Fisher) Platinum Hot Start PCR Master Mix (2X), 200 reactions
- 16S customized primer plates (see below)
- 1.7 ml tubes
- 0.5 ml Qubit tubes
- 96-well PCR plates
- Filter tips
- Loading dye
- Gel red dye
- SeaKem Agarose
- Hi-Lo DNA Marker

### Earth Microbiome Procedure:

This protocol is based on the 16S Earth Microbiome Protocol, which is available on the Earth Microbiome website: <http://www.earthmicrobiome.org>

#### 16S Amplification Primers

1. Use this protocol for the amplification of prokaryotes using paired-end 16S sequencing with primers 515f-806r, which target the V4 region of the 16S SSU rRNA (see Caporaso et al. 2012, *ISME*; Walters et al. 2016, *mSystems*).
2. Modified 515f-806r primer pair:
  - a. Barcodes are on the forward primer of the 515f-806rB primer pair.
  - b. Degeneracy was added to forward and reverse primers to remove biases (see Parada et al. 2014, *Environ Microbiol*; Apprill et al. 2015, *Aquat Microb Ecol*).

#### Primer Order

1. Order primers in the 5' -> 3' orientation from IDT.
2. Barcoded primers are ordered in 96-well plates at 100 nmol.
3. Check the option to receive same amount per well (8 nmol).

### Resuspension of Primers

1. Resuspend primers with ultra-pure water under sterile conditions in a hood.
2. Stock plate are 100 uM and should be aliquoted and diluted to 10 uM.
3. Use the standard desalting cleanup option.
4. If using a robot, plate volume should be 15 uL per well or greater.

### 515fB Forward Primer

1. 5' Illumina adapter (index sequencing primer)
2. Golay barcode
3. Forward primer pad
4. Forward primer linker
5. Forward primer (515fB)

AATGATACGGCGACCACCGAGATCTACACGCT XXXXXXXXXXXXX  
TATGGTAATT GT GTGYCAGCMGCCGCGGTAA

### 806rB Reverse Primer

1. Reverse complement of 3' Illumina adapter
2. Reverse primer pad
3. Reverse primer linker
4. Reverse primer (806rB)

CAAGCAGAAGACGGCATAACGAGAT AGTCAGCCAG CC  
GGACTACNVGGGTWTCTAAT

### 16S Amplification in Triplicate

1. Retrieve genomic material and put in fridge to thaw.
2. Get reagents for master mix from freezer.
3. Use 16S V4 515f-806rB primers in resuspended and aliquoted 96-well plates.
4. Let everything thaw on bench. Put primers back on ice after they have thawed.
5. Label 96-well plates. Samples are run in triplicate. Make plate maps of where each sample is placed in 96-well PCR plate.
6. Vortex & spin down everything. Put all components back on ice.
7. Set up PCR reaction:

PCR Reaction	Concentration	Volume
ddH <sub>2</sub> O		13 ul
16S V4 515f	0.2 uM	0.5 ul
16S V4 806rB	0.2 uM	0.5 ul
Platinum HotStart PCR Master Mix	0.8x	5 ul
DNA template		1 ul
Total		25 ul

8. Briefly centrifuge PCR plates.
9. Place PCR plates in Thermal Cycler. PCR should run for about two hours.

Cycle	1	35			1	1
Temp	94.0	94.0	50.0	72.0	72.0	4.0
Time	3:00	0:45	1:00	1:30	10:00	∞

### Pool PCR Triplicates

1. Pool PCR triplicates of each sample into a single volume (75 ul)
2. Create pool in fresh 96-well plate.

### Run Gel

1. Prepare 1.2% agarose gel and gel dye red.
2. Protocol: 40 minutes at 100V with DNA Hi-Lo Ladder.
3. Should yield 300-350 bp PCR product.
4. Bioanalyzer is an optional alternative method.

### Qubit & Pool Samples

1. Use Qubit protocol to record new concentrations of the samples.
2. Based on Qubit concentrations, combine an equimolar amount (240 ng) of amplicon from each sample in a single tube.
3. If spurious bands were present on gel, one-half of the final pool can be run on a gel and then gel extracted to select only the target bands.

### Clean Samples

1. Clean the amplicon pool following the instructions from Qiagen UltraClean PCR Cleanup Kit.
2. Qubit the cleaned sample before evaporation/dilution.

### MiSeq

1. Sample must be 2.0 nM to run on the MiSeq platform.
2. Use a V2 300 cycle kit.
3. 20% PhiX to increase diversity.